

# SURGISCAN RANGE OVERVIEW & SPECIFICATIONS

DataMatrix Readers	SurgiScan Ultra Lite	SurgiScan Ultra
<b>Basic Features</b>		
2D DataMatrix Reading	✓	✓
1D & QR Code Reading	✓	✗
Safe LED Lighting	✓	✓
Read Speed	1 read/sec	200 reads/sec
Minimum Code Size	0.8mm x 0.8mm	0.4mm x 0.4mm
Reading Distance Range	0-10cm	0-3cm
<b>Physical Characteristics</b>		
Weight	640g	640g
Reading Window Diameter	95mm	95mm
Casing Material	Aluminium	Aluminium
Sapphire Crystal Window	✓	✓
Identification Band Options	Black	✗
<b>Software Compatibility</b>		
Microsoft Windows	10 or 11	7, 10 or 11
Connection	USB 2.0 x1	USB 2.0 x1
<b>Data Transport Options</b>		
HTTP	✗	✓
Keyboard Wedge	✓	✓
Named Pipe	✗	✓
Serial (COM) Port	✓	✓
<b>Compatible Marking Methods</b>		
Laser Marked	✓	✓
Dot Peen	✓	✓
Chemical Etched	✓	✓
Microdot	✓	✓
Label	✓	✓
<b>Software Features</b>		
Live Video Window	✗	✓
Advanced Code Manipulation	✗	✓
Light Customisation	✗	✓
Code Quality Indicator	✗	✓
<b>Other Features</b>		
Benchtop	✓	✓
Handheld	✗	✗
Protective Caps Available	✓	✓
Plug & Play	✓	✗
<b>Recommended Retail Price</b>		
R.R.P	£995	£1,950

## Contact us for more information:

+44 (0) 1691 831 497

+44 (0) 1527 833 038

info@2dsurgical.com

3 Navigation Court  
Waterside  
Stoke Prior  
Bromsgrove  
United Kingdom  
B60 4FD